

Fluorescent Quantitative Detection system

LineGene K Series

User's Manual

Attention

Users are recommended to read the contents of this manual thoroughly before operating the Bioer Fluorescent Quantitative PCR Detection System.

To carefully observe all special Warnings and Cautions outlined in this manual.

This manual should be maintained properly in good condition for reference.

Caution: Copyright reserved. The Bioer Co. reserves the right to modify this manual at any time without notice.

The manual contains copyright protected and patented material. Without prior written consent from Hangzhou Bioer Technology Co., Ltd., any part of the manual shall not be duplicated, reproduced or translated into any other language.

***Thank you for your purchase of this product.
Before initial use of this instrument, please read this manual thoroughly !***

File No.: BYQ609600000ESM
File Version: Dec 2014 Version 1.0

Important Notes

1. Usual practice

Note: Very important information is contained within this manual and it should be carefully read before first use of the instrument. Failure to operate instrument according to the instruction could result in damage or abnormal functioning of the instrument.

Warning! The warning message requires extremely careful operation of a certain step. Failure to observe the instruction could result in serious personal injury.

2. Safety

During operation, maintenance and repair of this instrument, the following basic safety notes must be observed. In case of failure to follow these measures or the warnings or notes indicated herein, the basic protection provided by the instrument, its safety criteria of design and manufacture, and its predicted use range would be impaired.

Hangzhou Bioer Technology Co., Ltd. shall not be held responsible for any consequences resulting from the user's failure to observe the following requirements.

Note: The instrument, complying with the Standard GB4793.1/IEC61010-1, is a general instrument of class I, the protection degree is IP20. It is intended for indoor use of Elevation 200 meters below.

The instrument, complying with the Standard YY0648/IEC61010-2-101, is used for IVD medical equipment

a) Instrument earth

In order to avoid an electric shock, the input power cable of the instrument must be properly earthed. This instrument uses a 10A 3-core earthed plug, which is provided with a third (earth) pin. It is for use with an earth type power socket and is a safety unit. If the plug cannot be inserted into the socket, the socket must be fixed by a qualified electrician, to maintain the safety function of the plug and the protection it provides.

b) Keeping apart from the live circuit

Internal maintenance or replacement of any part of the thermal cyclers must only be carried out by qualified personnel. The instrument must be disconnected from the mains circuit prior to any maintenance being carried out.

c) Use of power supply

Before connecting to the mains and switching the instrument on, make sure the voltage is consistent with the instruments requirements (100-240V, 50/60Hz). The rated load for the power socket must not be less than the instruments maximum load of 600W

d) Power wire

The instrument is supplied with a power cable which should be used at all times when operating the instrument. If the power cable is damaged it should be replaced with a new one of the same specification. Care must be taken that the power cable does not get compressed or tightly bent and that it does not lie across areas where it may cause a trip hazard to personnel.

e) Insertion and withdrawal of power cable

At insertion and withdrawal of power cable, the back of the plug shall be firmly held with the hand. The plug must be completely and tightly inserted into the socket and must not be

removed by pulling the cable. The back of the plug should be grasped in the hand and pulled directly backwards to remove from the socket.

f) Placement of instrument

This instrument should not be positioned in a place where it is difficult to cut off the power supply.

This instrument should be placed in a low relative humidity (RH) and low dust environment well away from any water (e.g. sinks and pipes). The room should be well ventilated, and free from corrosive gas, or interference from a strong magnetic field. The instrument should not be placed in a wet or dusty location, but should be positioned on a sturdy, level and secure table appropriate to its weight.

The openings on this instrument are for ventilation purposes and in order to avoid overheating of the instrument they shall not be blocked or covered. When a single set or several sets of instruments are used, the space between its ventilation openings and the nearest object should not be less than 30cm.

Excessive environmental temperature would impair the test performance and could result in failure of the instrument. This instrument should not be used in locations subjected to direct sunlight or strong radiation or light source, as this could impair the fluorescence detection.

The instrument should be kept away from hot gas, furnaces, stoves and all other sources of heat.

When switched off, the power should also be switched off. If the instrument is not going to be used for a long time, the power should be switched off, the power plug withdrawn and the instrument covered with soft cloth or plastic film to prevent dust or foreign bodies entering the machine.

g) Notes during operation

During test, care shall be taken to prevent liquid from dropping onto the instrument. The castoff used in test, such as consumables, reagent, and so on, should be treated as require, and should not be thrown away or poured.

During test, if there is hazardous substances, user must be trained before using.

Hazardous substances, which has been used, should be coped with and saved according to derrection for use.

User, who operates the instrument, must be trained and has relevant quantification.

Caution: If any of the following should occur, you should immediately switch off the power supply, withdraw the power plug from the power socket, and contact the supplier to effect a repair: Repairs can only be carried out by suitably qualified engineers.

- Liquid gets inside the instrument.
- The instrument is rained upon or water is spilled over it.
- The instrument works abnormally, or generates an abnormal sound/s or generates a strange odour.
- The instrument is dropped or its casing is damaged.
- There is an obvious change in the function of the instrument.

Caution: When you deal with potential contagious matter such as flesh sample or reagent, which is likely to touch skin, protecting glove or other protecting measures are need to be used.

h) Transportation again

If transporting the instrument again, the instrument and its detection wells need be cleared wholly before transportation, and be disinfected by UV light. After clearing, the fixing pin shall be tightened to transportation well when the instrument stops running and power off so that motion component of the instrument interior can be fixed and motion component is avoided to be damaged during transportation. Introduction of the location from transportation well is seen in chapter of fixing pin backout.

Manual handling steps,

a) **Think before lifting/handling.** Plan the lift. Can handling aids be used? Where is the load going to be placed? Will help be needed with the load? Remove obstructions such as discarded wrapping materials. For a long lift, consider resting the load midway on a table or bench to change grip.

b) **Adopt a stable position.** The feet should be apart with one leg slightly forward to maintain balance (alongside the load, if it is on the ground). The worker should be prepared to move their feet during the lift to maintain their stability. Avoid tight clothing or unsuitable footwear, which may make this difficult.

c) **Get a good hold.** Where possible, the load should be hugged as close as possible to the body. This may be better than gripping it tightly with hands only. **Start in a good posture.** At the start of the lift, slight bending of the back, hips and knees is preferable to fully flexing the back (stooping) or fully flexing the hips and knees (squatting). **Don't flex the back any further while lifting.** This can happen if the legs begin to straighten before starting to raise the load.

d) **Keep the load close to the waist.** Keep the load close to the body for as long as possible while lifting. Keep the heaviest side of the load next to the body. If a close approach to the load is not possible, try to slide it towards the body before attempting to lift it.





e) **Avoid twisting the back or leaning sideways,** especially while the back is bent. Shoulders should be kept level and facing in the same direction as the hips. Turning by moving the feet is better than twisting and lifting at the same time.

f) **Keep the head up when handling.** Look ahead, not down at the load, once it has been held securely. **Move smoothly.** The load should not be jerked or snatched as this can make it harder to keep control and can increase the risk of injury. **Don't lift or handle more than can be easily managed.** There is a difference between what people can lift and what they can safely lift. If in doubt, seek advice or get help.

g) **Put down, then adjust.** If precise positioning of the load is necessary, put it down first, then slide it into the desired position.

i) Warning Sign

Warning identification

DANGER!		Place pasted this mark in instrument, is danger, if the instrument is used irrelevancy.
SCALDING!		Place pasted this mark in instrument, causes high temperature and is scalding during use.
BIOHAZARD		Place pasted this mark in instrument, caused biohazard during use.
PROTECTIVE EARTH		Protective earth is near to the place pasted this mark in instrument





Warning mark




Warning! This indicates a surface which will be hot during and immediately following the running of a programme. Contact with this metal area will cause burns.

Warning! During the use of the instrument the user may come into contact with biologically hazardous materials. Rules for safe handling of such materials must be followed. The operator must be appropriately trained.

j) Signs on the external packaging

Up		Correct position of transport package is vertical upwards.
Breakables		Be treated with care when transportation, there is breakables in transport package,
Fear rain		Transport package is not rained.
The limit of stacking layer		Maximum stacking layer of the same package is 2.

Temperature limit		Temperature limit that transport package should keep is form -20C to 55C.
-------------------	---	---

3. EMC Consideration

Note: This is a class A equipment, only suitable for use in establishments other than domestic, and those directly connected to a low voltage power supply network.

4. Maintenance of instrument

Any stains on the instrument can be cleaned with soft cloth soaked with a gentle cleaning solution.

Heat conductive oil medium should not be used in the block wells of this instrument.

Module should not be left open for any period of time as this may allow dust to enter the instrument.

Warning!

- When cleaning the instrument, the power should be turned off.
- The instrument surface should not be cleaned with corrosive cleaning agents.
- The instrument module includes precise optics, dust, foreign matter and residue should be avoided.

5. After-sales services

The warranty content and scope are shown in the warranty sheet.

Note:

- After unpacking, immediately check the goods against the packing list. If any parts are damaged or missing, please contact the supplier immediately.
- After qualification of acceptance, complete the product acceptance sheet and send (or fax) the copied sheet to the supplier for filing and maintenance.
- Before first use of the product, the user shall complete the instrument registration form and send to Hangzhou Bioer Technology Co., Ltd. to obtain the correct operation password.
- After unpacking, the packing box and packing materials should all be kept in case it is required for transportation or service in the future.
- In the event that a repair is required, the instrument must be disinfected before being sent to the repair department. A decontamination sheet should be completed and sent together with the instrument. These are available on request from your local supplier.
- It is recommended that service personnel disinfect the instrument on receipt in the service department, before commencing any scheduled work.
- Hangzhou Bioer Technology Co., Ltd. shall bear no liability in the event of any damage to the instrument occurring during transportation to the service department due to improper packaging.

Contents

CHAPTER 1 GENERAL DESCRIPTION	1
1. APPLICATIONS	1
2. FEATURES	1
3. NOTES TO MODEL AND DESCRIPTION	2
4. PERFORMANCE PARAMETERS	2
5. GENERAL DESCRIPTION OF FUNCTIONS OF THE SOFTWARE	3
CHAPTER 2 PREPARATIONS	4
1. TRANSPORT AND STORAGE CONDITIONS OF THE INSTRUMENT	4
2. NORMAL WORKING CONDITIONS	4
3. PREPARATIONS BEFORE STARTING THE INSTRUMENT	4
3.1 Connection of power wire and communication wire	4
4. SYSTEM INSTALLATION AND UNLOADING	4
4.1 System Installation	4
4.2 System Unloading	5
CHAPTER 3 START	6
1. CHECKS BEFORE START	6
2. START	6
3. STARTING SOFTWARE INTERFACE	7
CHAPTER 4 ABSOLUTE QUANTIFICATION	8
1. DESIGN EXPERIMENT	8
1.1 Create New Absolute Quantitative Experiment	8
1.2 Detector Setting	9
1.3 Sample Information Setting	10
1.4 Reaction Plate Setting	11
1.5 Programme Setting	13
2. PREPARE FOR REACTION	14
3. RUN THE EXPERIMENT	15
3.1 Preparation for reagent sample	15
3.2 Run Fluorescence Curve	16
3.3 Run Temperature Curve	17
3.4 Programme Setting	18
3.5 Working state indication lamps on instrument	18
3.5.1. Prompts which may occur during running	19
4. EXPERIMENT ANALYSIS	20
4.1 Check Results	20
4.1.1 Check the Amplification Plot	20
4.1.2 Check Standard Curve	23
4.1.3 Check Melting Curve	24
4.2 Adjusting Parameters and Re-analysis	25
5. EXPERIMENT REPORT	26
5.1 Designing a Report Template	26
5.2 Print Setting	27
5.3 Comprehensive Report	27
5.4 Report Printing	28
5.5 QC Summary	28
6. DATA EXPORT	29
6.1 Export to Database	29
6.2 Experiment Filing	29

6.3 Export Experiment Data to EXCEL	29
6.4 Export Experiment Data to TEXT	29
CHAPTER 5 RELATIVE QUANTITATIVE	30
1. DESIGN EXPERIMENT	30
1.1 Create New Relative Quantitative Experiment	30
1.2 Detector Setting	30
1.3 Sample Information Setting	31
1.4 Reaction Plate Setting	33
1.5 Programme Setting	34
2. PREPARE FOR REACTION	36
3. RUN THE EXPERIMENT	37
3.1 Run Fluorescence Curve	37
3.2 Run Temperature Curve	39
3.3 Programme Setting	40
4. EXPERIMENT ANALYSIS	41
4.1 Check Results	41
4.1.1 Check the Amplification Plot	41
4.1.2 Check Standard Curve	42
4.2 Check Relative Quantification	43
4.3 Adjust Parameter Reanalysis	44
5. EXPERIMENT REPORT	46
5.1 Comprehensive Report	46
5.2 QC Summary	47
6. DATA EXPORT	48
6.1 Export to Database	48
6.2 Experiment Filing	48
6.3 Export Experiment Data to EXCEL	48
6.4 Export Experiment Data to TEXT	48
CHAPTER 6 SNP	49
1. DESIGN EXPERIMENT	49
1.1 Create SNP Experiment	49
1.2 Detector Setting	49
1.3 Sample Information Setting	51
1.4 Reaction Plate Setting	52
1.5 Programme Setting	53
2. PREPARE FOR REACTION	55
3. RUN THE EXPERIMENT	56
3.1. Run Fluorescence Curve	56
3.2 Run Temperature Curve	58
3.3 Programme Setting	59
4. EXPERIMENT ANALYSIS	60
4.1 Check Results	60
4.1.1 Check the Amplification Plot	60
4.1.2 Check SNP	62
4.2 Adjust Parameter Re-analysis	63
5. EXPERIMENT REPORT	64
5.1 Designing a Report Template	64
5.2 Print Setting	65
5.3 Comprehensive Report	65
5.4 Report Printing	66
5.5 QC Summary	66
6. DATA EXPORT	67
6.1 Export to Database	67
6.2 Experiment Filing	67

6.3 Export Experiment Data to EXCEL	67
6.4 Export Experiment Data to TEXT	67
CHAPTER 7 HIGH RESOLUTION MELTING	68
1. DESIGN EXPERIMENT	68
1.1 Create High Resolution Melting Experiment	68
1.2 Detector Setting	68
1.3 Sample Information Setting	70
1.4 Reaction Plate Setting	71
1.5 Programme Setting	72
2. PREPARE FOR REACTION	74
3. RUN THE EXPERIMENT	75
3.1. Run Fluorescence Curve	75
3.2 Run Temperature Curve	76
3.3 Programme Setting	78
4. EXPERIMENT ANALYSIS	79
4.1 Check Results	79
4.1.1 Check the Amplification Plot	79
4.1.2 Check the Standard Curve	81
4.1.3 Check HRM	82
4.2 Adjust Parameter Re-analysis	84
5. EXPERIMENT REPORT	85
5.1 Comprehensive Report	85
5.2 QC Summary	86
6. DATA EXPORT	87
6.1 Export to Database	87
6.2 Experiment Filing	87
6.3 Export Experiment Data to EXCEL	87
6.4 Export Experiment Data to TEXT	87
CHAPTER 8 SERVICE	88
1. USER MANAGEMENT	88
2. EXPERIMENT MANAGEMENT	89
2.1 Experiment Management	89
2.2 Deleted Experiment Management	90
3. TEMPLATE MANAGEMENT	91
4. USER LOGIN	91
5. CHANGE PASSWORD	91
6. SEE RUNNING EXPERIMENT	92
CHAPTER 9 TOOL USAGE	93
1. GAIN SETTING	93
2. BLOCK SCAN METHOD	93
3. DETECTOR LIBRARY	93
4. CUSTOMIZED DYES	94
5. CUSTOMIZE COLUMNS	94
6. COLUMN SELECTION	95
7. SAMPLE COLUMN LIBRARY	96
8. INSTRUMENT CALIBRATION PARAMETERS	96
9. MEASURE CROSSTALK CALIBRATION PARAMETERS	97
10. CROSSTALK GAIN PARAMETER MEASUREMENT	97
11. SYSTEM MAINTENANCE	98
12. UPGRADE EXPERIMENT FILE FORMAT	98
13. TA CALCULATOR	99
CHAPTER 10 OTHER FUNCTIONS	100

1. INSTRUMENT OPERATION	100
1.1 Connect.....	100
1.2 Disconnect	100
1.3 Instrument Information.....	100
2. DATA QUERY	100
3. SYSTEM HELP.....	101
CHAPTER 11 MAINTENANCE	102
1. REGULAR CLEANING.....	102
2. ANALYSIS AND TROUBLESHOOTING	102
APPENDIX: LINEGENE K SERIES WIRING	105

Chapter 1 General description

This chapter mainly describes the applications, features, specification, model, performance parameters and software functions of this 48-well Fluorescent Quantitative Detection System.

1. Applications

This detection system adopts fluorescent real-time detection method to analyse PCR template amplification and is suitable for polymerase chain reaction fluorescent quantitative detection in clinical diagnosis fields of virus, tumor, hereditary disease.


2. Features

- Novel and human-orientated running interface for smooth operation.
- The adopted fluorescent real-time detection mode realizes simultaneous amplification and detection in the same tube without need of post experimental treatment.
- Advanced thermoelectric technology ensures fast and steady heating and cooling of the ultra-fast heat cycling system.
- Multi-point temperature control ensures consistent temperature of 48 sample wells.
- It can create a temperature gradient with 4 Thermo electric modules.
- The constant temperature function of SOAK enables low-temperature storage of PCR reagents.
- It uses maintenance-free long life LED excitation light source.
- The advanced fibre optic transmission technology makes the photo-electric detection system very sensitive and reliable.
- Precise optical path system and ultra-sensitive PMT system provide the most accurate and sensitive fluorescent detection.
- It can create a real-time dynamic monitor of the entire PCR amplification process.
- It has high linear range up to 10 orders of start DNA copies without serial dilution.
- It is unnecessary to open the PCR reaction tube, ensuring samples are protected from contamination during and after PCR and ensures accurate results.
- Multiplexing is possible.
- The hot-lid technology allows for oil-free operation of PCR.
- Automatic hot-lid technology needs no manual opening/closing and ensures constant pressure of the hot-lid used with different height reaction tubes or plates.
- User friendly interface with flexible programme setting and analysis and reporting using the stored parameters.
- It can print out one or more sample report(s).
- Automatic, correct and timely remote networked services provide the latest technical support.
- Unique patented advanced module bottom fluorescent detection technology.

3. Notes to model and description


Model:

FQD – 48 A



Description:

LineGene K plus



4. Performance parameters

Specification/model	FQD-48A							
Sample capacity	48x0.2ml (suitable for single tube, 8-row tubes and 48-well fully-skirted plate)							
Detection channel	F1	F2	F3	F4				
Applicable dye	FAM, SYBR Green I	VIC, HEX, TET, JOE, Cy3,	ROX, TEXAS -RED, TAMR A	Cy5 Quasar -670				
Temperature range of block working	4~105°C (Minimum division:0.1°C)							
Heating/cooling rate	4.0°C/s (max)							
Temperature fluctuation	≤±0.1°C (full-range), (55°C typical value ≤±0.1°C)							
Temperature accuracy	≤±0.2°C (full-range), (55°C typical value ≤±0.1°C)							
Temperature uniformity	≤±0.4°C (full-range), (55°C typical value ≤±0.3°C)							
Temperature range of hot-lid working	70~110°C (adjustable, default 105°C)							
Repeatability of fluorescent intensity detection	5%							
Running mode	Continuous running							
Operation system	Windows XP/Windows Vista/Windows7/Windows8							
Power supply	100-240V~ 50/60Hz 600W							
Dimensions	384mmx353mmx348 mm							
Weight	12kg							

5. General description of functions of the software

a) Parameters setting-up function (including temperature, time, cycles, heating/cooling rate, selection of detection channel and yield of photo-electric amplification tube).

b) Note function of text contents.

c) Sample material record function (sample No., sample name and sample data).

d) Document running display function (PCR heat cycle data display, fluorescence detection data display and real-time display of each data during running of instrument).

e) Detection data analysis function (The analysis function may be independently used without connection to the instrument).

f) Analysis result output function. It may output the analysis result to various types of document, e.g.: EXCEL, TXT document. It is possible to run an enquiry and print out analysis result, modify the printing format and select/de-select items to print.

g) Document storage function (setting up data, running data and analysis results).

h) Fault protection and alarm function

Caution: The above-mentioned software functions are merely for reference. The software functions may be modified without notice.

Chapter 2 Preparations

This chapter mainly describes use, transport and storage condition, structural composition, installation/unloading of software and preparations before first using the LineGene K series fluorescent quantitative detection system.

1. Transport and storage conditions of the instrument

Environmental temperature: $-20^{\circ}\text{C} \sim 55^{\circ}\text{C}$

Relative Humidity: $\leq 80\%$

2. Normal working conditions

Environmental temperature: $10^{\circ}\text{C} \sim 30^{\circ}\text{C}$

Environmental RH: $\leq 70\%$

Altitude: $< 2000\text{ m}$

Pollution degree: 2

Power supply: $100\text{-}240\text{V} \sim 50/60\text{Hz}$ 600W

Caution: Before using the instrument, please make sure the working conditions meets the above requirements. The power socket shall be a 3-hole socket and properly earthed.

3. Preparations before starting the instrument

3.1 Connection of power wire and communication wire

Connection of power wire: Only the power cable supplied with the instrument should be used. At connection, make sure the instrument power switch is in the "OFF" position and after connection the power cable should be checked to ensure a tight contact with the instrument socket; otherwise it should be replaced.

Connection of communication wire: The supplied communication wire and communication conversion box should be used. One end of the communication conversion box is connected with the DB15 communication interface at the back of the instrument, and the other end connected with the computer USB, RS232C or Bluetooth interface. After connection, screws should be securely tightened.

Caution: If with prolonged use the supplied power cable connection becomes loose, it should be replaced with one of the same type and specification.

The communication conversion box is built with special circuits and must not be opened.

4. System Installation and Unloading

4.1 System Installation

System Environment

Operating system: Windows XP/Windows Vista/Windows7/Windows8

Runtime environment: Net Framework 4.0

Other software: PDF reader

Minimum configuration:

Processor: Intel Core i3

Memory: 2GB

Hard Disc: 10GB

System installation

Double click PcrServer installation file (PcrServerSetup.exe) ► display the installation interface (select installation language) ► set up installation path ► install

Double click LineGene K plus installation file (LineGene K plus Setup.exe) ► display the installation interface (select installation language) ► set up installation

Operating system

1. Double click the LineGene K shortcut on the desktop Or Click the start menu ► Program ► LineGene K plus
2. Double click the PcrServer shortcut on the desktop Or Click the start menu ► Program ► PcrServer

4.2 System Unloading

Control Panel ► Add/Delete Program ► PcrServer ► Unloading

Control Panel ► Add/Delete Program ► LineGene K Plus ► Unloading

Chapter 3 Start

1. Checks before start

After inserting the power plug and switching on this detection system, the following should be checked:

- Check the voltage of the power supply is consistent with the system-required one.
- Check the power cable plug for correct and reliable insertion into the power socket.
- Check the communication converter for correct insertion into the host and tightening and the cable plug for correct and reliable connection to the computer.
- Check the update shift switch MODE of communication converter is set to “normal” state.
- Check the environmental conditions meet the required tolerances.

2. Start

In order to ensure effective connection and communication between the instrument and the computer system, the system shall be started in the following sequence:

1st step: Start computer display and host.


2nd step: Turn on power using the switch on the back of the instrument.

3rd step: Press the Run Switch on the front cover of the instrument to set the system ready to run.

4th step: On the computer after entering operation system, start the LineGene K fluorescent quantitative detection system. To start the software, click “LineGene K plus” from [Start]/[programme] menu or double click the short-cut icon on desktop.

3. Starting software interface


Run the programme and it will display the system window.



The system window consists of the menu bar, the toolbar and the main page.

Chapter 4 Absolute Quantification

1. Design Experiment




This section describes how to design a new absolute quantification experiment and covers inspection item setting, sample information setting, reaction plate setting and programme setting.

1.1 Create New Absolute Quantitative Experiment

1. Click build **Absolute** on the **Home** interface and this will open the absolute quantitative experiment window.


NOTE: The Absolute quantitative experiment can be also created by:

- Clicking **File ► New ► Absolute** on the menu bar
- Clicking **New ► Absolute** on the toolbar



1.2 Detector Setting

1. Click **Setup** ► **Detector**



2. Input experiment properties

Input the experiment name, user name and any comments in the experiment properties column.



A screenshot of the 'Experiment Properties' dialog box. It contains two text input fields: 'Experiment Name' with the value '20111117_Experiment' and 'User Name' with the value 'user'. To the right, there is a 'Comment' field with the value 'remark'.


3. Detector Setting

Set up the Detector, Assay, Dye and Colour.



If necessary, the user can also:

- a. Add detector
- b. Add assay
- c. Delete detector
- d. Delete assay
- e. Add the detector in the Detector Library: click **Add Detector From Library** ► the **Detector Library** window will pop up ► select the Detector in the window to be added

The user can also conduct Add, Modify and Delete operations in the item library.



- f. Set up the detector, set up the assay, set up the dye name and set up the colour

Detector	Reporter	Color	Master Mix	Primer	Probe	Supplies	Batch Number
Target1	FAM						
Target2	FAM						

4. Set up reference dye

Reference Dye

VIC

1.3 Sample Information Setting

1. Click **Setup** ► **Sample**

Setup

- Detector
- Sample
- Plate
- Program

2. Add sample information

a. Itemized addition: input ID in **Sample ID** ► press **Enter** ► add information for one sample

b. Batch addition: click **Batch Add** ► the Batch Add window will pop up

Batch Add

Start Sample Id Sample Count

Add Cancel

3. Delete sample information

a. Itemized deletion: select one sample ► click **Delete** ► delete the selected sample information

b. Delete all: click **Clear All** ► delete all sample information






4. Import/Export sample information

a. Click **Import Sample Info** ► the File Import window will pop up ► import sample information file in CSV format

b. Click **Export Sample Info** ► the Save As window will pop up ► the sample information will be exported in CSV file format


Sample ID Batch Add Delete Clear All Import Samples Info Export Samples Info

5. Set up sample information

Samples					
Sample Id	Color	Sample Name	Sampling Time	Submitting Date	
a1		Sample1	2013-12-06	2013-12-06	
a2		Sample2	2013-12-06	2013-12-06	
a3		Sample3	2013-12-06	2013-12-06	
a4		Sample4	2013-12-06	2013-12-06	
a5		Sample5	2013-12-06	2013-12-06	

1.4 Reaction Plate Setting


1. Click **Setup** ► **Plate**







2. Set up the inspection criteria of the reaction plate

- Select reaction plate well site: click Reaction Plate well Site

The user can also right click the reaction plate well site to Copy, Paste and Add New Detector. Adding a new detector will open the **Edit Detector Library** window.



b. Select Assay item and modify the property, concentration and concentration unit.


Property	Name	Concentration	Concentration unit
	Unknown	NO	Copies/ml
	Standard	YES	IU/ml
	Negative	NO	Fg/ml
	Positive	NO	Pg/ml

c. Select a sample and the list displayed will change

d. Zoom-In, Zoom-Out and reset the reaction plate.

e. Sample Auto Arrange


f. Check Well Table



The screenshot displays a software interface for assay configuration and well table management. On the left, the 'Detectors' panel shows 'Assay Item' set to 'Target1 - FAM', 'Property' set to 'U', and 'Con.' checked. Below it, the 'Concentration Unit' is set to 'copies/ml'. The 'Samples' panel shows a list of samples (A1-A5) with 'A1' selected. On the right, the 'Well Table' is a grid with columns for '#', 'Well', 'Sample ID', 'Sample Name', 'Assay Item', 'Property', 'Ct SD', and 'Con. SD'. The first row is populated with '1', 'A01', 'A1', and 'Target1'. At the bottom, a reaction plate grid shows wells A1 through A24, with well A1 highlighted and containing a 'U' icon.

1.5 Programme Setting

1. Click Setup ► Programme



2. Run Programme Setup

a. Create new stage: the user can create a new **Hold Stage**, **Cycling Stage** or **Melting Stage**

The user can also click **Add Stage** directly and the default will be creating a new **Cycling Stage**.

b. Create new step: the user can create a new step **Before** or **After** the currently selected step


The user can also click **Add Step** and the default will be adding a new Step at the end of the currently selected stage or after the currently selected step.

c. Delete: the user can delete the currently selected step or stage


d. Display form: click **Display With Table** ► new window will pop up ► the details of the current experiment will be displayed in a table.

e. Set up the experimental data of the hold stage, cycling stage and melting stage melting section

f. Set up the hot-lid temperature and liquid volume




2. Prepare for Reaction



The user should make full preparations prior to the experiment:

- Ensure appropriate materials are used.
- Ensure the arrangement of the PCR reaction plate is consistent with the setting layout of the reaction plate in Section 1.4.

3. Run the Experiment



This section describes how to run/operate the experiment after loading the reaction plate and includes how to operate the fluorescence curve, the temperature curve and programming

Caution: Before starting, make checks and follow the procedure for correct start up of the system. The green lamp of the run switch will be lit and the system will be ready to run.


Caution: Before running the programme, push the module smoothly until the locking sound is heard and the alarm lamp is switched off. If the module is pulled out or is improperly closed, the software will produce a pop-up warning and the alarm lamp will light up. In this case, the temperature programme can be run but the fluorescence scanning data would become invalid.

3.1 Preparation for reagent sample

- Preparation for reagent: The LineGene K series fluorescent quantitative PCR detection system uses 0.2ml PCR tubes, strip tubes or 48 well PCR plates to conduct the reaction. The recommended reaction volume is 10µl~50µl for an optimal reaction system.
- The tube, strip tube or 48 well PCR plate, must have an optically clear bottom.
- Centrifugal operation: Before placing reactions into the instrument, it is recommended that a short centrifugal spin is used to ensure that the reagent is at the bottom of the reaction tube and the reagent/sample mix is free from bubbles.
- To insert test tube: If individual tubes or strip tubes are used, and the sample quantity is less than the maximum capacity of the instrument, it is recommended the sample tubes should be evenly distributed across the block as far as possible. This will create even pressure across the hot-lid, ensuring consistent pressure on all the tubes during running, which greatly improves temperature consistency across all sample tubes.



Correct.
The sample is at the bottom of the PCR tube




Incorrect


1. Requires a greater spin speed
2. Requires a longer spin time

3.2 Run Fluorescence Curve

1. Click Run ► Fluorescence Curve




2. Click Start Run



3. Operating confirmation


- Modify hot-lid temperature and liquid quantity (sample volume).
- Gain (baseline) parameter setting
- Target fluorescence value setting



4. After it starts operating, the user can:


- Skip the current stage
- Add a cycle
- Delete a cycle
- Stop run

5. Plot display setting
 - a. Assay item
 - b. Plot colour




3.3 Run Temperature Curve

1. Click Run ► Temperature Curve




2. Click Start Run




3. Operating confirmation

- a. Modify hot-lid temperature and liquid quantity (sample volume).
- b. Gain (baseline) parameter setting
- c. Target fluorescence value setting



4. After it starts running, the user can:
- Skip the current stage
 - Add a cycle
 - Delete a cycle
 - Stop run



3.4 Programme Setting

The user can only check the programme setting but cannot make modifications.

3.5 Working state indication lamps on instrument

The panel at the right of the instrument is fixed with 1 lamp and the colors related to the system state during the running of a programme:

- **Standby:** The indicator lamp lights **blue**, which denotes that the entire machine is ready to operate.
- **Running:** The indicator lamp lights **green**, which denotes that the entire machine is running a programme.
- **Error:** The indicator lamp lights **red**, which denotes that the instrument has detected a fault.

Caution: For prolonged shutdown, switch off the power at the back of the instrument and at the socket. When switched on again, the hot-lid and module will revert to the default settings.

The front cover of the instrument is fixed with a self-locking key to control energizing of its internal control system:

- Run Switch: running/standby switch.
- After pressing this key, the green indicator lamp is lit on the instrument, the internal system is energized and the instrument is ready to run the programme.
- After pressing this key again, the key will spring out, the green indicator lamp goes off, the instrument internal system is de-energized and the system is under standby.


Caution: The run switch is for ease of operation and is merely used for temporary or short term closing down of the control system. When the system is under the standby state, the instrument internal AC circuit remains live.

3.5.1. Prompts which may occur during running

- Hot-lid temperature sensor alarm prompt
- Sink temperature sensor alarm prompt
- Environmental temperature sensor alarm prompt
- Module temperature sensor alarm prompt
- Module sensor short-circuit or short-circuit alarm prompt

Caution: In case the temperature alarm displays during the running of a programme, the PCR detection system will terminate the current programme. The instrument should be switched off and then re-started.

4. Experiment Analysis




This section describes how to view the experiment analysis results after running an experiment and adjusting parameters for re-analysis.

This section covers the analysis of amplification curves and standard curves, adjusting parameters for re-analysis and importing parameters.

4.1 Check Results


4.1.1 Check the Amplification Plot

1. Click **Analysis** ► **Amplification Plot**



2. Check the amplification curve
 - a. Set up colour
 - b. Set up plot type
 - c. Set up show dye

When the background colour of a dye name is blue, it will be displayed; while white indicates it will not be displayed.



3. Check the reaction plate

- a. Select reaction plate well site and check corresponding well site curve
The default is all wells are selected
- b. Zoom-In, Zoom-Out and reset the reaction plate
- c. Check well table
- d. Check results summary




4. Set up assay

- a. Set up assay
- b. Set up threshold
- c. Set up automatic baseline


When the threshold value is not automatic, the user cannot set up the automatic

baseline




4.1.2 Check Standard Curve

1. Click Analysis ► Standard Curve



2. Check standard curve

a. Set up assay



3. Check the reaction plate
 - a. Select reaction plate well site and check corresponding well site curve
The default is all wells are selected
 - b. Zoom-In, Zoom-Out and reset the reaction plate
 - c. Check well table information
 - d. Check results summary




4.1.3 Check Melting Curve

1. Click **Analysis** ► **Melting Curve**



2. Check the melting curve
 - a. Check the fluorescence curve
 - b. Check the derivative curve
 - c. Set up colour



3. Check the reaction plate
 - a. Select reaction plate well site and check corresponding well site curve
The default is all wells are selected
 - b. Zoom-In, Zoom-Out and reset the reaction plate
 - c. Check well table information
 - d. Check results summary




4. Set up assay
 - a. Set up assay
 - b. Set up colour




4.2 Adjusting Parameters and Re-analysis

1. Click **Analysis Settings** ► the Analysis Settings dialog box will pop up
 - a. Adjust the start cycle and end cycle of the baseline
 - b. Adjust Ct analysis algorithm
 - c. Set up the use of S fitting
 - d. Set up the stage to use for Ct analysis
 - e. Set up the automatic threshold value
 - f. Advanced setting
 - g. Standard curve setting



5. Experiment Report



This section describes how to print an experiment report and covers designing of a report template and print settings.

5.1 Designing a Report Template

1. Click **Report** ► **Report Template Editor** ► the Report Designer window will pop up

The report consists of controls and the user can add, modify and delete controls. Available controls include Static Text, Dynamic Text, Line, Static Image, Amplification curve and Quantification Analysis Results.

The screenshot shows the Report Designer window with the following components:

- Available controls:**
 - Common Controls: Static Text, Dynamic Text, Static Image, Line, Amplification Curve, Quantification Analysis Result
 - Known Controls: Static Text Controls, Dynamic Text Controls
- Appearance:**
 - Alignment: MiddleRight
 - BackColor: White
 - Border: Solid, 1, False, False, False
 - Color: Black
 - Font: Tahoma, 8.25pt
- Data:** Tag
- Design:** DesignVisible: True, Name: Label10
- Layout:** Location: 93, 62, Padding: 0, 0, 0, 0, Size: 100, 20, Type: Label

The report template layout includes:


- Header: [Hospital]
- Sub-header: [Report]
- Form fields: Name: [Name], Sex: [Sex], Age: [Age], Hospital No.: [Hospital No.]
- Table:

Test Item	Test Result	Reference	Conclusion
- Graph: Amplification Curve (Y-axis: 0 to 5000, X-axis: 0 to 45)
- Footer: [Submitting Date], Report Date: [ReportDate], Tester: [Tester], Checker: [Checker]

5.2 Print Setting

1. Click **Report** ► **Print Template Setting** ► the Print Template Setting window will open

The user can set up the laboratory name, report name, reference value, tester, checker, amplification plot, default report template and paper size.



The dialog box is titled "Print Template Settings (Absolute)". It contains several sections:


- Template Setup:** Fields for Hospital, Report, Reference (set to 100), Tester, and Checker.
- Amplification Plot Setup:** Legend options: Color, LineStyle.
- Print Setup:** Default Report Template (dropdown, set to default) and Paper Size (dropdown, set to A4).
- Printer:** Use Default Printer, Use Custom Printer (with a dropdown menu).

Buttons for "OK" and "Cancel" are at the bottom right.

5.3 Comprehensive Report

1. Click **Report** ► **Consolidated Reports** ► the Consolidated Report window will pop up

The Consolidated Report includes the basic information, sample information, amplification curve, standard curve, plate information, etc..



The "Consolidated Report" window displays the following information:

Experiment Name: 20141030_131639
Experiment Type: Absolute
User Name: admin
File Name: C:\Bioer\LineGene4800 for research\4840-英文\experiment\扩增20141030_131639.fqd
Run Time: 2014/10/30 13:17:47 - 2014/10/30 14:44:58
Gain: F1:7, F2:19, F3:19, F4:7

Run Program

Hold Stage

Target	Incubation Time	Rate	Sampling
94	60	4	<input type="checkbox"/>

PCR Stage Cycles:45

Target	Incubation Time	Rate	2nd Temp.	Step Size	Step Delay	Grad Temp.	Grad Range	Sampling
94	10	4						<input type="checkbox"/>
55	30	4						<input checked="" type="checkbox"/>

Detectors

Detector	Reporter	Color	Master Mix	Primer	Probe	Supplies	Batch Number
Target1	FAM	Yellow					
	HEX	Blue					
	ROX	Green					
	CY5	Magenta					

Plot Plate

Table Plate

#	Well	Assay Item	Property	Dye	Std. Con.	Sample ID
1	A01	Target1	Unknown	FAM		


Report Items

- Basic Information
- Run Program
- Detectors
- Plot Plate
- Table Plate
- Amplification Curve
- Plot Type: Linear
- Quan. Analysis Result

Buttons: Create Report, Print Report


5.4 Report Printing

1. Click Report ► Report Print




2. Report print setting

- Set up report template
- Print setting (please refer to Section 5.2)
- Select items to print
- Print preview
- Print the report




5.5 QC Summary


1. Click Report ► QC Summary



2. Check the QC summary



6. Data Export



This section describes how to export data and covers exporting to a database, experiment filing and exporting the experiment data to EXCEL.

6.1 Export to Database

Click **Data Summary** ► **Export to Database** ► the Save File dialog box will pop up ► save the exported database file

6.2 Experiment Filing

1. Click **Data Summary** ► **Archived Experiment Directory** ► the Experimental archive storage directory window will pop up ► set up the storage path of file.



2. Experiment filing

Click **Data Summary** ► **Archived Experiment** ► export the filed experiment file
The suffix of the filed experiment file is .fqh

6.3 Export Experiment Data to EXCEL


Click **Data Summary** ► **Export Experiment** ► **Export Experiment to Excle** ► the exported experiment data will generate EXCEL file

6.4 Export Experiment Data to TEXT

Click **Data Summary** ► **Export Experiment** ► **Export Experiment to Text** ► the exported experiment data will generate TEXT file

Chapter 5 Relative Quantitative

1. Design Experiment




This section describes how to design a relative quantitative experiment and covers creating new relative quantitative experiment, inspection item setting, sample information setting, reaction plate setting and programme setting.

1.1 Create New Relative Quantitative Experiment


1. Click **Relative** on **Home** interface and create Relative Quantitative Experiment window. Relative quantitative experiment can be also created by:

- Clicking **New ► Relative** on the toolbar
- Clicking **File ► New ► Relative** on the menu bar



1.2 Detector Setting

1. Click **Setup ► Detector**



2. Input Experiment Properties

Input the Experiment name, User name and Comment in the basic information column.

3. Inspection Item Setting



a. Set up the Detector, Assay, Dye and Colour.

b. Add detector

c. Delete detector

d. Add detector from library

The user can also conduct Add, Modify and Delete operations in the item library.



4. Set up reference dye



1.3 Sample Information Setting

1. Click **Setup** ► **Sample**



2. Add sample information

- a. Itemized addition: input ID in **Sample ID** ► press **Enter** ► add information for one sample.
- b. Batch addition: click **Batch Add** ► the Batch Add window will pop up




3. Delete sample information

- a. Itemized deletion: select one sample ► click **Delete** ► delete the selected sample information
- b. Delete all: click **Clear All** ► delete all sample information

4. Import/Export sample information

- a. Click **Import Sample Info** ► the File Import window will pop up ► import sample information file in CSV format
- b. Click **Export Sample Info** ► the Save As window will pop up ► the sample information will be exported in CSV file format



5. Set up sample information

Samples				
Sample Id	Color	Sample Name	Sampling Time	Submitting Date
a1			2013-12-06	2013-12-06
a2			2013-12-06	2013-12-06
a3			2013-12-06	2013-12-06
a4			2013-12-06	2013-12-06
a5			2013-12-06	2013-12-06

1.4 Reaction Plate Setting


1. Click **Setup** ► **Plate**






2. Set up the inspection criteria of the reaction plate

a. Select reaction plate well site: click Reaction Plate well Site

The user can also right click the reaction plate well site to Copy, Paste and Add New Detector. Adding a new detector will open the **Edit Detector Library** window.



b. Select inspection item and modify the property, concentration and concentration unit.

Property	Name	Concentration	Concentration unit
	Unknown	NO	Copies/ml
	Standard	YES	IU/ml
	Negative	NO	Fg/ml
			Pg/ml

c. Select a sample and the list displayed will change

d. Zoom-In, Zoom-Out and reset the reaction plate.

e. Sample Auto Arrange

f. Check Well Table

1.5 Programme Setting

1. Click **Setup** ► **Programme**



2. Run Programme Setup

- a. Create new stage: the user can create a new **Hold Stage**, **Cycling Stage** or **Melting Stage**

The user can also click **Add Stage** directly and the default will be creating a new **Cycling Stage**.

- b. Create new step: the user can create a new step **Before** or **After** the currently selected step


The user can also click **Add Step** and the default will be adding a new Step at the end of the currently selected stage or after the currently selected step.

- c. Delete: the user can delete the currently selected step or stage


- d. Display form: click **Display With Table** ► new window will pop up ► the details of the current experiment will be displayed in a table.

e. Set up the experimental data of the hold stage, cycling stage and melting stage
melting section

f. Set up the hot-lid temperature and liquid volume




2. Prepare for Reaction



The user should make full preparations prior to the experiment

- Ensure appropriate materials are used.
- Ensure the arrangement of the PCR reaction plate is consistent with the setting layout of the reaction plate in Section 1.4.


3. Run the Experiment




This section describes how to run/operate the experiment after loading the reaction plate and includes how to operate the fluorescence curve, the temperature curve and programming

3.1 Run Fluorescence Curve

1. Click Run ► Fluorescence Curve




2. Click Start Run



3. Operating confirmation

- Modify hot-lid temperature and liquid quantity (sample volume)
- Gain (baseline) parameter setting
- Target fluorescence value setting



4. After it starts running, the user can:

- a. Skip the current stage
- b. Add a cycle
- c. Delete a cycle
- d. Stop run

5. Plot display setting

- a. Assay item
- b. Plot colour

Experiment Type: **Relative** Experiment Name: **20150127_084452**

Run Status

Fluorescence Curve

Fluorescence

Sampling Count

Plot Settings


Assay Item: ALL Color: Well

Well Table


#	Well	Sample ID	Sample Name	Assay Item	Property	Ct	SD
1	A01			Target1	Unknown		
1	A01			Target2	Unknown		
2	A02						
3	A03						
4	A04			Target1	Unknown		
4	A04			Target2	Unknown		
5	A05						
6	A06						
7	A07			Target1	Standard		
7	A07			Target2	Standard		
8	A08						
9	A09						
10	A10			Target1	Unknown		
10	A10			Target2	Unknown		
11	A11						
12	A12						

3.2 Run Temperature Curve

1. Click Run ► Temperature Curve




2. Click Run ► Start




3. Operating confirmation

- Modify hot-lid temperature and liquid quantity (sample volume)
- Gain (baseline) parameter setting
- Target fluorescence value setting



4. After it starts running, the user can:


- Skip the current stage
- Add a cycle
- Delete a cycle
- Stop run



3.3 Programme Setting

The user can only check the programme setting but cannot make modifications.

4. Experiment Analysis



This section describes how to view the experiment analysis results after running an experiment and adjusting parameters for re-analysis.

This section covers the analysis of amplification curves and standard curves, the analysis of relative quantification, adjusting parameters for re-analysis and importing parameters.

4.1 Check Results

4.1.1 Check the Amplification Plot


1. Click **Analysis** ► **Amplification Plot**



2. Check the amplification curve

- a. Set up colour
- b. Set up plot type
- c. Set up show dye

When the background colour of a dye name is blue, it will be displayed; while white indicates it will not be displayed.



3. Check the reaction plate

- a. Select reaction plate well site and check corresponding well site curve
The default is all wells are selected
- b. Zoom-In, Zoom-Out and reset the reaction plate
- c. Check well table
- d. Check results summary




4. Set up assay

- a. Set up assay
- b. Set up threshold
- c. Set up automatic baseline


When the threshold value is not automatic, the user cannot set up the automatic

Baseline




4.1.2 Check Standard Curve

1. Click Analysis ► Standard Curve



2. Check standard curve
 a. Set up assay




3. Check the reaction plate

- Select reaction plate well site and check corresponding well site curve
 The default is all wells are selected
- Zoom-In, Zoom-Out and reset the reaction plate
- Check well table
- Check results summary




4.2 Check Relative Quantification

1. Click Analysis ► Relative Quantification



2. Check relative quantitative

- Set up the show type



b. Check the analysis results

Sample Id	Assay Item	Property	GOI Aver. Con.	GOI Con. SD	HKG Aver. Con.	HKG Con. SD	Max	Min	Aver.
01	target1	Comparison	7.99e+03	0.00e+00	1.37e+04	0.00e+00	1	1	1
02	target1	Unknown	1.10e+07	1.05e+06	1.93e+05	1.48e+04	63.92	49.95	56.94
03	target1	Unknown	8.48e+05	1.31e+05	6.14e+04	9.61e+03	16.84	10.78	13.81
04	target1	Unknown	9.40e+04	1.40e+04	3.67e+04	2.06e+03	2.97	2.15	2.56
06	target1	Unknown	3.72e+03	2.66e+01	3.08e+04	8.82e+02	0.12	0.12	0.12
07	target1	Unknown	9.44e+05	1.43e+05	3.95e+04	6.33e+03	29.18	18.63	23.9
08	target1	Unknown	9.33e+04	3.53e+04	2.73e+04	5.86e+03	4.9	1.93	3.41
09	target1	Unknown	4.14e+03	2.62e+03	2.33e+04	8.42e+02	0.29	0.07	0.18
11	target1	Unknown	8.44e+06	5.34e+05	9.71e+04	3.93e+04	122.5	51.28	86.89
12	target1	Unknown	7.21e+04	1.20e+03	1.57e+04	2.97e+02	4.7	4.47	4.58
13	target1	Unknown	1.10e+04	0.00e+00	1.51e+04	0.00e+00	0.73	0.73	0.73
14	target1	Unknown	8.12e+06	8.33e+05	8.05e+04	1.74e+04	125.02	76.77	100.89
16	target1	Unknown	8.25e+05	6.25e+04	2.50e+04	2.87e+03	37.59	28.5	33.05
16	target1	Unknown	6.87e+03	3.28e+03	8.01e+03	4.28e+02	1.27	0.45	0.86

4.3 Adjust Parameter Reanalysis

1. Click **Analysis Settings** ► the Analysis Settings dialog box will pop up
 - a. Adjust the start cycle and end cycle of the baseline
 - b. Adjust Ct analysis algorithm
 - c. Set up the use of S fitting
 - d. Set up the stage to use for Ct analysis
 - e. Set up the automatic threshold value
 - f. Advanced setting
 - g. Standard curve setting
 - h. Relative quantification setting

Analysis Settings

Ct Settings **Advanced Settings** Standard Curve Settings Relative Quantification Settings

The stage to use for Ct analysis:

The algorithm to calculate Ct: S Fitting

Assay Item	Threshold	Start Cycle	End Cycle
target1 - SYBR	230.6	Auto	Auto
target1 - HEX	Auto	Auto	Auto

target1 - SYBR

Auto Threshold


Threshold:

Auto Baseline

Start Cycle: End Cycle:

5. Experiment Report

This section describes how to print experiment report and covers report template designing and print setting.



This section describes how to print an experiment report and covers designing of a report template and print settings.

5.1 Comprehensive Report

1. Click **Report** ► **Consolidated Reports** ► the Consolidated Report window will pop up
The Consolidated Report includes the basic information, sample information, amplification curve, standard curve, plate information, etc..

Consolidated Report 1 / 5

Experiment Name: 20141124_080143
 Experiment Type: Relative
 User Name: admin
 File Name: I:\48PLUS\400101\相对定量_20141124_080143.fqd
 Run Time: 2014/11/24 08:04:21 - 2014/11/24 09:28:54
 Gain: F1:12, F2:14
 Reference Sample:

Run Program

Target	Incubation Time	Rate	Sampling
94	60	4	<input type="checkbox"/>

PCR Stage Cycles:45

Target	Incubation Time	Rate	2nd Temp.	Step Size	Step Delay	Grad Temp.	Grad Range	Sampling
94	10	4						<input type="checkbox"/>
55	30	4						<input checked="" type="checkbox"/>

Detectors

Detector	Endogenous Control	Reporter	Color	Master Mix	Primer	Probe	Supplies	Batch Number
项目1		FAM	GREEN					
项目2		HEX	MAGENTA					

Plot Plate

#	Well	Assay Item	Property	Dye	Std. Con.	Sample ID
1	A01	项目1	Unknown	FAM		

Report Items


- Basic Information
- Run Program
- Detectors
- Plot Plate
- Table Plate
- Amplification Curve
 - Plot Type: Linear
- Quan. Analysis Result
- Relative Plot
 - Y Axis: Linear
 - X Axis: Sample
- Relative Analysis Result

Create Report


Print Report

5.2 QC Summary


1. Click Report ► QC Summary



2. Check the QC summary



6. Data Export




This section describes how to export data and covers exporting to a database, experiment filing and exporting the experiment data to EXCEL.

6.1 Export to Database

Click **Data Summary** ► **Export to Database** ► the Save File dialog box will pop up ► save the exported database file

6.2 Experiment Filing

1. Click **Data Summary** ► **Archived Experiment Directory** ► the Experimental archive storage directory window will pop up ► set up the storage path of file



2. Experiment filing

Click **Data Summary** ► **Archived Experiment** ► export the filed experiment file
The suffix of the filed experiment file is .fqh

6.3 Export Experiment Data to EXCEL

Click **Data Summary** ► **Export Experiment** ► **Export Experiment to Excle** ► the exported experiment data will generate EXCEL file

6.4 Export Experiment Data to TEXT

Click **Data Summary** ► **Export Experiment** ► **Export Experiment to Text** ► the exported experiment data will generate TEXT file

Chapter 6 SNP

1. Design Experiment

Start
▼
Design experiment
▼
Prepare for the reaction
▼
Run the experiment
▼
Experiment analysis
▼
Experiment report
▼
Data export
▼
End


This section describes how to design an SNP experiment and covers creating a new SNP experiment, inspection item setting, sample information setting, reaction plate setting and programme setting.

1.1 Create SNP Experiment

1. Click **SNP** on **Home** interface and create SNP Experiment window.


An SNP experiment can be also created by:

- a. Clicking **New ▶ SNP** on the toolbar
- b. Clicking **File ▶ New ▶ SNP** on the menu bar



1.2 Detector Setting

1. Click **Setup ▶ Detector**



2. Input basic information

Input the experiment name, user name and any comments in the experiment properties column.


3. Inspection Item Setting

Set up the Detector, Allele, Dye and Colour.


If necessary, the user can also:

- a. Add Detector
- b. Delete Detector
- c. Add the Detector in the Detector library: click **Add Detector From Library** ► the **Detector Library** window will pop up ► select the Detector on the window to be added

The user can also conduct Add, Modify and Delete operations in the item library.



- d. Set up the item name, set up the dye name and set up the colour




4. Set up reference dye



1.3 Sample Information Setting


1. Click **Setup** ► **Sample**



2. Add sample information

a. Itemized addition: input ID in **Sample ID** ► press **Enter** ► add information for one sample

b. Batch addition: click **Batch Add** ► the Batch Add window will pop up



3. Delete sample information


a. Itemized deletion: select one sample ► click **Delete** ► delete the selected sample information

b. Delete all: click **Clear All** ► delete all sample information

4. Import/Export sample information

a. Click **Import Sample Info** ► the File Import window will pop up ► import sample information file in CSV format

b. Click **Export Sample Info** ► the Save As window will pop up ► the sample information will be exported in CSV file format




5. Set up sample information

Sample Id	Color	Sample Name	Sampling Time	Submitting Date	
a1		Sample1	2013-12-06	2013-12-06	
a2		Sample2	2013-12-06	2013-12-06	
a3		Sample3	2013-12-06	2013-12-06	
a4		Sample4	2013-12-06	2013-12-06	
a5		Sample5	2013-12-06	2013-12-06	

1.4 Reaction Plate Setting


1. Click Setup ► Plate








2. Set up the inspection criteria of the reaction plate

a. Select reaction plate well site: click Reaction Plate well Site

The user can also right click the reaction plate well site to Copy, Paste and Add New Detector. Adding a new detector will open the **Edit Detector Library** window.



b. Select inspection item and modify the property, concentration and concentration unit.

Property	Name	Concentration	Concentration unit
	Unknown	NO	Copies/ml IU/ml Fg/ml Pg/ml
	Negative	NO	
	Positive Allelic gene 1	NO	
	Positive Heterozygous	NO	
	Positive Allelic gene 2	NO	

c. Select a sample and the list displayed will change

d. Zoom-In, Zoom-Out and reset the reaction plate.


e. Sample Auto Arrange

f. Check Well Table

The screenshot displays the software's configuration and data management interface. On the left, the 'Detectors' panel is set to 'Target1' with a 'U' property and 'copies/ml' unit. Below it, the 'Samples' panel lists five samples (a1-a5) with color-coded icons. The 'Well Table' on the right is a table with columns: #, Well, Sample ID, Sample Name, Assay Item, Property, and Con. SD. The first two rows are populated with '1 A01 a1 Target1 Unknown'. At the bottom, a 24-well plate grid is shown with wells A1 and B1 highlighted.

1.5 Programme Setting

1. Click Setup ► Programme



2. Run Programme Setup

- a. Create new stage: the user can create a new **Hold Stage**, **Cycling Stage** or **Melting Stage**

The user can also click **Add Stage** directly and the default will be creating a new **Cycling Stage**.

- b. Create new step: the user can create a new step **Before** or **After** the currently selected step


The user can also click **Add Step** and the default will be adding a new Step at the end of the currently selected stage or after the currently selected step.

- c. Delete: the user can delete the currently selected step or stage


d. Display form: click **Display With Table** ► new window will pop up ► the details of the current experiment will be displayed in a table.

e. Set up the experimental data of the hold stage, cycling stage and melting stage melting section

f. Set up the hot-lid temperature and liquid volume



2. Prepare for Reaction




The user should make full preparations prior to the experiment:

- Ensure appropriate materials are used.
- Ensure the arrangement of the PCR reaction plate is consistent with the setting layout of reaction plate in Section 1.4.

3. Run the Experiment


This section describes how to run the experiment after loading the reaction plate and covers the operating of fluorescence curve, the operating of temperature curve and programme setting.




This section describes how to run/operate the experiment after loading the reaction plate and includes how to operate the fluorescence curve, the temperature curve and programming

3.1. Run Fluorescence Curve

1. Click Run ► Fluorescence Curve




2. Click Start Run



3. Operating confirmation

- Modify hot-lid temperature and liquid quantity (sample volume).
- Gain (baseline) parameter setting
- Target fluorescence value setting



4. After it starts running, the user can:
 - a. Skip the current stage
 - b. Add a cycle
 - c. Delete a cycle
 - d. Stop run
5. Plot display setting
 - a. Assay item
 - b. Plot colour

Experiment Type: **SNP** Experiment Name: **20150127_091219**

Run Status

Fluorescence Curve


Plot Settings

Assay Item: Target1 Color: Well


#	Well	Sample ID	Sample Name	Assay Item	Property
1	A01				
2	A02				
3	A03			Target1	Positive Allele
4	A03			Target1	Positive Allele
5	A04			Target1	Positive Allele
6	A04			Target1	Positive Allele
7	A05			Target1	Positive Allele
8	A05			Target1	Positive Allele
9	A06			Target1	Positive Allele
10	A06			Target1	Positive Allele
11	A07				
12	A07				
13	A08				
14	A08				
15	A09				
16	A09				
17	A10				
18	A10				
19	A11			Target1	Positive Allele
20	A11			Target1	Positive Allele

3.2 Run Temperature Curve

1. Click Run ► Temperature Curve




2. Click Start Run




3. Operating confirmation

- Modify hot-lid temperature and liquid quantity (sample volume).
- Gain (baseline) parameter setting
- Target fluorescence value setting



4. After it starts running, the user can:


- Skip the current stage
- Add a cycle
- Delete a cycle
- Stop run



3.3 Programme Setting

The user can only check the programme setting but cannot make modifications.

4. Experiment Analysis




This section describes how to view the experiment analysis results after running an experiment and adjusting parameters for re-analysis.

This section covers the analysis of amplification curves and standard curves, adjusting parameters for re-analysis and importing parameters.

4.1 Check Results

4.1.1 Check the Amplification Plot


1. Click **Analysis** ► **Amplification Plot**



2. Check the amplification curve

- a. Set up colour
- b. Set up plot type
- c. Set up show dye

When the background colour of a dye name is blue, it will be displayed; while white indicates it will not be displayed.




3. Check the reaction plate
 - a. Select reaction plate well site and check corresponding well site curve
The default is all wells are selected
 - b. Zoom-In, Zoom-Out and reset the reaction plate
 - c. Check well table
 - d. Check results summary



4. Set up inspection item

- a. Set up assay
- b. Set up threshold
- c. Set up automatic baseline

When the threshold value is not automatic, the user cannot set up the automatic baseline



4.1.2 Check SNP

1. Click Analysis ► SNP




2. Check SNP

- a. Select well site

The user can select well site by dragging a rectangle with the mouse around the wells of interest or select wells one by one.


- b. Set up Assay
- c. Set up manual calls




3. Check the reaction plate
 - a. Select reaction plate well site and check corresponding well site curve
The default is all wells are selected
 - b. Zoom-In, Zoom-Out and reset the reaction plate
 - c. Check well table information
 - d. Check results summary

4.2 Adjust Parameter Re-analysis

1. Click **Analysis Settings** ► the Analysis Settings dialog box will pop up
 - a. Adjust analysis data
 - b. Adjust whether the inspection item will retain manual recognition genotype



5. Experiment Report




This section describes how to print an experiment report and covers designing of a report template and print setting.

5.1 Designing a Report Template

1. Click **Report** ► **Report Template Editor** ► the Report Designer window will pop up


The report consists of controls and the user can add, modify and delete controls. Available controls include Static Text, Dynamic Text, Line, Static Image and SNP Typing Curve.



5.2 Print Setting

1. Click Report ► Print Template Setting ► the Print Template Setting window will pop up

The user can set up the laboratory name, report name, reference value, tester, checker, amplification plot set up, default report template and paper size.



The dialog box is titled "Print Template Settings (SNP)" and contains the following sections:


- Template Setup:** Includes input fields for Hospital, Report, Tester, and Checker.
- Print Setup:** Includes a dropdown for Default Report Template (set to "default") and a dropdown for Paper Size (set to "A4").
- Printer:** Includes radio buttons for "Use Default Printer" (selected) and "Use Custom Printer" (with a dropdown menu).

Buttons for "OK" and "Cancel" are located at the bottom right.

5.3 Comprehensive Report

1. Click Report ► Consolidated Reports ► the Consolidated Report window will pop up.

The Consolidated Report includes the basic information, sample information, amplification curve, SNP, plate information, etc..



The "Consolidated Report" window displays the following information:

Experiment Information:
 Experiment Name: 20141218_110546
 Experiment Type: SNP
 User Name: admin
 File Name: I:\48PLUS\400101\SNP .fqd
 Run Time: 2014/12/18 11:06:47 - 2014/12/18 11:18:11
 Gain: F1:7, F2:7

Run Program:
 PCR Stage Cycles:40

Target	Incubation Time	Rate	2nd Temp.	Step Size	Step Delay	Grad Temp.	Grad Range	Sampling
60	10	4						<input checked="" type="checkbox"/>

Detectors:

Detector	Allele	Reporter	Color	Master Mix	Primer	Probe	Supplies	Batch Number
项目1	Allele1	FAM	Green					
	Allele2	HEX	Red					

Plot Plate: A grid showing amplification curves for 24 wells (A-D, 1-6). Wells 11-14 and 21-24 show positive results with FAM (green) and HEX (red) peaks.


Table Plate:

#	Well	Assay Item	Property	Dye	Std. Con.	Sample ID
3	A03	项目1	Positive Allele 1	FAM		
3	A03	项目1	Positive Allele 1	HEX		
4	A04	项目1	Positive Allele 1	FAM		
4	A04	项目1	Positive Allele 1	HEX		
5	A05	项目1	Positive Allele 1	FAM		
5	A05	项目1	Positive Allele 1	HEX		

Report Items: A sidebar on the right lists report components with checkboxes: Basic Information, Run Program, Detectors, Plot Plate, Table Plate, Amplification Curve, Quant. Analysis Result, SNP Typing Plot, and SNP Analysis Result. The "Plot Type" is set to "Linear". Buttons for "Create Report" and "Print Report" are at the bottom.


5.4 Report Printing

1. Click Report ► Report Print




2. Report print setting

- Set up report template
- Print setting (please refer to Section 5.2)
- Select print items
- Print preview
- Print the report




5.5 QC Summary


1. Click Report ► QC Summary



2. Check the QC summary



6. Data Export



This section describes how to export data and covers exporting to a database, experiment filing and exporting the experiment data to EXCEL.

6.1 Export to Database

Click **Data Summary** ► **Export to Database** ► the Save File dialog box will pop up ► save the exported database file

6.2 Experiment Filing

1. Click **Data Summary** ► **Archived Experiment Directory** ► the Experimental archive storage directory window will pop up ► set up the storage path of file



2. Experiment filing

Click **Data Summary** ► **Archived Experiment** ► export the filed experiment file

The suffix of filed experiment file is .fqh

6.3 Export Experiment Data to EXCEL


Click **Data Summary** ► **Export Experiment** ► **Export Experiment to Excle** ► the exported experiment data will generate EXCEL file.

6.4 Export Experiment Data to TEXT

Click **Data Summary** ► **Export Experiment** ► **Export Experiment to Text** ► the exported experiment data will generate TEXT file.

Chapter 7 High Resolution Melting

1. Design Experiment




This section describes how to design an SNP experiment and covers creating a new SNP experiment, inspection item setting, sample information setting, reaction plate setting and programme setting.

1.1 Create High Resolution Melting Experiment

1. Click **HRM** on **Home** interface and create SNP Experiment window.


An SNP experiment can be also created by:

- Clicking **New ► HRM** on the toolbar
- Clicking **File ► New ► HRM** on the menu bar



1.2 Detector Setting

1. Click **Setup ► Detector**



2. Input basic information

Input the experiment name, user name and any comments in the experiment properties column.


3. Inspection Item Setting

Set up the Detector, Allele, Dye and Colour.


If necessary, the user can also:

- a. Add Detector
- b. Delete Detector
- c. Add the Detector in the Detector library: click **Add Detector From Library** ► the **Detector Library** window will pop up ► select the Detector on the window to be added

The user can also conduct Add, Modify and Delete operations in the item library.



- d. Set up the item name, set up the dye name and set up the colour




4. Set up reference dye



1.3 Sample Information Setting


1. Click **Setup** ► **Sample**



2. Add sample information

a. Itemized addition: input ID in **Sample ID** ► press **Enter** ► add information for one sample

b. Batch addition: click **Batch Add** ► the Batch Add window will pop up



3. Delete sample information


a. Itemized deletion: select one sample ► click **Delete** ► delete the selected sample information

b. Delete all: click **Clear All** ► delete all sample information

4. Import/Export sample information

a. Click **Import Sample Info** ► the File Import window will pop up ► import sample information file in CSV format

b. Click **Export Sample Info** ► the Save As window will pop up ► the sample information will be exported in CSV file format




5. Set up sample information

A screenshot of a "Samples" table with 6 columns: Sample Id, Color, Sample Name, Sampling Time, and Submitting Date. The table contains 5 rows of sample data.

Sample Id	Color	Sample Name	Sampling Time	Submitting Date
a1	Red	Sample1	2013-12-06	2013-12-06
a2	Green	Sample2	2013-12-06	2013-12-06
a3	Black	Sample3	2013-12-06	2013-12-06
a4	Blue	Sample4	2013-12-06	2013-12-06
a5	Purple	Sample5	2013-12-06	2013-12-06

1.4 Reaction Plate Setting


1. Click Setup ► Plate



2. Set up the inspection criteria of the reaction plate

- a. Select reaction plate well site: click Reaction Plate well Site

The user can also right click the reaction plate well site to Copy, Paste and Add New Detector. Adding a new detector will open the **Edit Detector Library** window.



- b. Select inspection item and modify the property, concentration and concentration unit.

Property	Name	Concentration	Concentration unit
U	Unknown	NO	Copies/ml
S	Standard	YES	IU/ml
N	Negative	NO	Fg/ml
P	Positive	NO	Pg/ml

- c. Select a sample and the list displayed will change
- d. Zoom-In, Zoom-Out and reset the reaction plate.
- e. Sample Auto Arrange
- f. Check Well Table

The screenshot displays a software interface for configuring a PCR assay. It is divided into three main sections:

- Detectors:** A table with columns 'Assay Item', 'Property', and 'Con.'. It lists 'Target1 - FAM' and 'Target1 - HEX'. Below this table is a 'Concentration Unit' dropdown set to 'copies/ml'.
- Samples:** A table with columns 'Sample ID' and 'Sample Name'. It lists samples 'a1' through 'a5', each with a colored icon and a checkbox.
- Well Table:** A large grid with columns: '#', 'Well', 'Sample ID', 'Sample Name', 'Assay Item', 'Property', 'Ct SD', and 'Con. SD'. The first two rows are populated with data for wells A01 and B01.

At the bottom of the interface is a 2x24 grid representing a microplate layout, with columns numbered 1-24 and rows labeled A and B. A gear icon is visible in the bottom right corner.

1.5 Programme Setting

1. Click **Setup** ► **Programme**

The 'Setup' menu is shown with the following options:

- Setup (dropdown arrow)
- Detector (right arrow)
- Sample (right arrow)
- Plate (right arrow)
- Program (right arrow)

2. Run Programme Setup

- a. Create new stage: the user can create a new **Hold Stage**, **Cycling Stage** or **Melting Stage**

The user can also click **Add Stage** directly and the default will be creating a new **Cycling Stage**.

- b. Create new step: the user can create a new step **Before** or **After** the currently selected step


The user can also click **Add Step** and the default will be adding a new Step at the end of the currently selected stage or after the currently selected step.

c. Delete: the user can delete the currently selected step or stage


d. Display form: click **Display With Table** ► new window will pop up ► the details of the current experiment will be displayed in a table.

e. Set up the experimental data of the hold stage, cycling stage and melting stage melting section

f. Set up the hot-lid temperature and liquid volume



2. Prepare for Reaction




The user should make full preparations prior to the experiment:

- Ensure appropriate materials are used.
- Ensure the arrangement of the PCR reaction plate is consistent with the setting layout of reaction plate in Section 1.4.

3. Run the Experiment


This section describes how to run the experiment after loading the reaction plate and covers the operating of fluorescence curve, the operating of temperature curve and programme setting.




This section describes how to run/operate the experiment after loading the reaction plate and includes how to operate the fluorescence curve, the temperature curve and programming

3.1. Run Fluorescence Curve

1. Click Run ► Fluorescence Curve




2. Click Start Run




3. Operating confirmation

- Modify hot-lid temperature and liquid quantity (sample volume).
- Gain (baseline) parameter setting
- Target fluorescence value setting




4. After it starts running, the user can:
 - a. Skip the current stage
 - b. Add a cycle
 - c. Delete a cycle
 - d. Stop run
5. Plot display setting
 - a. Assay item
 - b. Plot colour




3.2 Run Temperature Curve

1. Click Run ► Temperature Curve




2. Click **Start Run**




3. Operating confirmation

- a. Modify hot-lid temperature and liquid quantity (sample volume).
- b. Gain (baseline) parameter setting
- c. Target fluorescence value setting



4. After it starts running, the user can:


- a. Skip the current stage
- b. Add a cycle
- c. Delete a cycle
- d. Stop run



3.3 Programme Setting

The user can only check the programme setting but cannot make modifications.

4. Experiment Analysis




This section describes how to view the experiment analysis results after running an experiment and adjusting parameters for re-analysis.

This section covers the analysis of amplification curves and standard curves, adjusting parameters for re-analysis and importing parameters.

4.1 Check Results

4.1.1 Check the Amplification Plot


1. Click **Analysis** ► **Amplification Plot**



2. Check the amplification curve

- a. Set up colour
- b. Set up plot type
- c. Set up show dye

When the background colour of a dye name is blue, it will be displayed; while white indicates it will not be displayed.



3. Check the reaction plate
 - a. Select reaction plate well site and check corresponding well site curve
The default is all wells are selected
 - b. Zoom-In, Zoom-Out and reset the reaction plate
 - c. Check well table
 - d. Check results summary




4. Set up inspection item
 - a. Set up assay
 - b. Set up threshold
 - c. Set up automatic baseline
When the threshold value is not automatic, the user cannot set up the automatic baseline

Assay

Assay: Target1-FAM Threshold: Auto 267 Auto Baseline


4.1.2 Check the Standard Curve

1. Click Analysis ► Standard Curve



2. Check the Standard Curve

a. Set up array




3. Check the reaction plate

- Select reaction plate well site and check corresponding well site curve
The default is all wells are selected
- Zoom-In, Zoom-Out and reset the reaction plate
- Check well table
- Check results summary




4.1.3 Check HRM

1. Click Analysis ► HRM Curve




2. Check the fluorescence curve

- a. Set up target
- b. Set up color




3. Check the derivative curve


- a. Set up target
- b. Set up color



4. Check the aligned curve
 - a. Set up target
 - b. Set up color



5. Check the Different Pilot
 - a. Set up target
 - b. Set up color




6. Check the reaction plate
 - a. Select reaction plate well site and check corresponding well site curve
The default is all wells are selected
 - b. Zoom-In, Zoom-Out and reset the reaction plate
 - c. Check well table




4.2 Adjust Parameter Re-analysis

1. Click **Analysis Settings** ► the Analysis Settings dialog box will pop up
 - a. Adjust analysis data
 - b. Adjust whether the inspection item will retain manual recognition genotype



5. Experiment Report



This section describes how to print an experiment report and covers designing of a report template and print setting.

5.1 Comprehensive Report

1. Click **Report** ► **Consolidated Reports** ► the Consolidated Report window will pop up.

The Consolidated Report includes the basic information, sample information, amplification curve, HRM curve, plate information, etc..

Plot Plate

Well	Assay Item	Property	Dye	Std. Con.	Sample ID
A10	项目 1	Unknown	SYBR		
A11	项目 1	Unknown	SYBR		
A12	项目 1	Unknown	SYBR		
A18	项目 1	Unknown	SYBR		
A19	项目 1	Unknown	SYBR		
A20	项目 1	Unknown	SYBR		
B05	项目 1	Unknown	SYBR		
B06	项目 1	Unknown	SYBR		
B07	项目 1	Unknown	SYBR		
B14	项目 1	Unknown	SYBR		

Table Plate

#	Well	Assay Item	Property	Dye	Std. Con.	Sample ID
10	A10	项目 1	Unknown	SYBR		
11	A11	项目 1	Unknown	SYBR		
12	A12	项目 1	Unknown	SYBR		
18	A18	项目 1	Unknown	SYBR		
19	A19	项目 1	Unknown	SYBR		
20	A20	项目 1	Unknown	SYBR		
29	B05	项目 1	Unknown	SYBR		
30	B06	项目 1	Unknown	SYBR		
31	B07	项目 1	Unknown	SYBR		
38	B14	项目 1	Unknown	SYBR		

Report Items

- Basic Information
- Run Program
- Detectors
- Plot Plate
- Table Plate
- Melting Curve
- Melting Curve(Derivative)
- Melting Analysis Result

Create Report

Print Report

5.2 QC Summary

1. Click Report ► QC Summary

Report ▼

Report Print ▶

QC Summary ▶

2. Check the QC summary

Amplification Plot

Color Well Plot Type Linear

Fluorescence

Cycle

A01	A02	A03	A04	A05	A06	A07	A08	A09	A10	A11	A12	A13	A14	A15	A16	A17	A18	A19	A20	A21
B01	B02	B03	B04	B05	B06	B07	B08	B09	B10	B11	B12	B13	B14	B15	B16	B17	B18	B19	B20	B21

QC Summary

Description	Value	Use	Result
Negative control with a Ct less than	38	<input checked="" type="checkbox"/>	
Positive control with a Ct greater than	30	<input checked="" type="checkbox"/>	
Unknown without a Ct	N/A	<input checked="" type="checkbox"/>	
Standard without a Ct	N/A	<input checked="" type="checkbox"/>	

6. Data Export

Start
▼
Design experiment
▼
Prepare for the reaction
▼
Run the experiment
▼
Experiment analysis
▼
Experiment report
▼
Data export
▼
End

This section describes how to export data and covers exporting to a database, experiment filing and exporting the experiment data to EXCEL.

6.1 Export to Database

Click **Data Summary** ► **Export to Database** ► the Save File dialog box will pop up ► save the exported database file

6.2 Experiment Filing

1. Click **Data Summary** ► **Archived Experiment Directory** ► the Experimental archive storage directory window will pop up ► set up the storage path of file



2. Experiment filing

Click **Data Summary** ► **Archived Experiment** ► export the filed experiment file

The suffix of filed experiment file is .fqh

6.3 Export Experiment Data to EXCEL

Click **Data Summary** ► **Export Experiment** ► **Export Experiment to Excle** ► the exported experiment data will generate EXCEL file.

6.4 Export Experiment Data to TEXT


Click **Data Summary** ► **Export Experiment** ► **Export Experiment to Text** ► the exported experiment data will generate TEXT file.

Chapter 8 Service

1. User Management

User management is used to manage user information

Click **Service ► User Management** on the menu bar




The screenshot shows a window titled "User Management" with a table of users. The table has columns for User ID, User Name, Disabled, Locked, CreateTime, Last Login Time, Last Login Address, Last Locked Time, Failed Count, and Permission. There are three rows of data.

User ID	User Name	Disabled	Locked	CreateTime	Last Login Time	Last Login Address	Last Locked Time	Failed Count	Permission
1	admin	NO	NO	2014-04-29	2014-05-06	127.0.0.1	0001-01-01	0	Manage experiment, Run experiment, View experiment, Use
2	aaa	NO	NO	2014-05-04				0	Manage experiment, Run experiment, View experiment, Use
3	bbb	NO	NO	2014-05-04				0	


The user can:

- a. delete user
- b. update user




The screenshot shows a dialog box titled "Update user". It contains a text field for "User name:" with the value "aaa". Below it are radio buttons for "Disabled:" with "Not disabled account" selected. There are four checked checkboxes for "Permission:": "User manage", "Run experiment", "View experiment", and "Manage experiment". At the bottom are "OK" and "Cancel" buttons.

c. change password



A dialog box titled "Change password" with a close button (X) in the top right corner. It contains three text input fields: "Old password:", "New password:", and "Confirm new password:". Below the fields are three buttons: "OK", "Reset", and "Cancel".

d. add user



A dialog box titled "Add user" with a close button (X) in the top right corner. It contains three text input fields: "User name:", "Password:", and "Confirm password:". Below the fields is a "Permission:" section with four checkboxes: "User manage", "Run experiment", "View experiment", and "Manage experiment". At the bottom are two buttons: "OK" and "Cancel".


2. Experiment Management

Experiment Management is used to manage experiment information and deleted experiment information.

2.1 Experiment Management

Click **Service ► Experiment management ► Experiment management** on the menu bar the user can:

- a. clear query condition
- b. set query condition
- c. query
- d. delete experiment
- e. download experiment
- f. edit experiment




2.2 Deleted Experiment Management

Click **Service** ► **Experiment Management** ► **Deleted Experiment Management** on the menu bar

The user can:

- a. clear query condition
- b. set query condition
- c. query
- d. delete experiment
- e. recover experiment
- f. clear experiment




3. Template Management

Template Management is used to manage template information.

Click **Service ► Template Management** on the menu bar


The user can:

- a. download template
- b. delete template




4. User Login

Click **Service ► User Login** on the menu bar



5. Change Password

Click **Service ► Change Password** on the menu bar



6. See Running Experiment

See Running Experiment is used to see running experiment which is running on connected instrument.

Click **Service ► See Running Experiment** on the menu bar

The screenshot shows the software interface for "See Running Experiment". At the top, it displays "Experiment Type: Absolute" and "Experiment Name: 20150129_090152". The main area is divided into several sections:

- Run Status:** Contains a "Fluorescence Curve" plot with "Fluorescence" on the y-axis (0 to 1000) and "Sampling Count" on the x-axis (0 to 40).
- Plot Settings:** Includes dropdown menus for "Assay Item" (set to "Target1-FAM") and "Color" (set to "Well").
- Well Table:** A table with the following data:

#	Well	Sample ID	Sample Name	Assay Item	Property	Ct SD
1	A01	A1		Target1	Unknown	
2	A02					
3	A03					
4	A04					
5	A05					
6	A06					
7	A07					
8	A08					
9	A09					
10	A10					
11	A11					
12	A12					
13	A13					
14	A14					
15	A15					
16	A16					
- Well Plate Visualization:** A 2x24 grid of wells. The first well (A1) is highlighted in blue and contains a small icon.

Chapter 9 Tool Usage


1. Gain Setting

The **Gain Setting** tool is used to set up gain modes.

Click **Tools** ► **Gain Setting** ► the following window will pop up

Gain setting can be set up as: **reference gain, custom gain and auto gain**


In **Custom Gain** mode, the user can modify the gain value.



2. Block Scan Method

Click **Tools** ► **Block Scan Method** ► the following window will pop up.

The user can select Whole Block Scan or Line Scan.




3. Detector Library

The **Detector Library** tool is used to set up the inspection libraries of absolute quantitative, relative quantitative and SNP analysis.

Click **Tools** ► **Detector Library** ► **(Absolute /Relative/SNP)** ► open the following window

The user can:

- a. Add Detector
- b. Modify Detector
- c. Delete Detector



4. Customized Dyes


The **Customized Dyes** tool is used to set up existing dyes and newly added dyes.

Click **Tools** ► **Customize Dyes** ► open the following window

The user can:

- a. Create dye
- b. Modify dye name and channel
- c. Delete dye
- d. Move dye upward
- e. Move dye downward

After adding new dyes or modifying dyes, the user should conduct crosstalk parameter measurements.




5. Customize Columns

Click **Tools** ► **Customize Columns** ► the following window will pop up

The user can:

- a. Add columns
- b. Delete columns
- c. Modify column name




6. Column Selection

The **Select Columns** tool is used to add the new columns in above section into current existing columns, or remove existing columns in current column.

Click **Tools** ► **Select Columns** ► the following window will pop up

1. Current existing column items include sample, report, report setting, query and query condition
2. Double click column can add or remove a column
3. Column with (*) indicates it cannot be removed




7. Sample Column Library

The **Sample Column Library** tool is used in the experiment design phase. The user can select the definition of contents in the drop-down box when setting up sample information.

Click **Tools** ► **Sample Column Library** ► the following window will pop up

The user can:


- Add columns
- Delete columns
- Edit the columns content



8. Instrument Calibration Parameters

The **Instrument Calibration Parameters** tool is used to calibrate the instrument parameters.

Click **Tools** ► **Instrument Calibration Parameters** ► the following window will pop up



9. Measure Crosstalk Calibration Parameters

The **Measure Crosstalk Calibration Parameters** tool is used to measure crosstalk correction parameters.

Click **Tools ► Measure Crosstalk Calibration Parameters ►** the following window will pop up

The user can add and modify the channels to be tested and dyes according to his needs; upload corresponding reaction plates and operate the experiment. When the experiment is over, the system will automatically save the crosstalk correction parameters.

Experiment Type: **Absolute** Experiment Name: **Crosstalk Parameter Measurement**

Experiment Properties

Experiment Name: Comment:

User Name:

Detectors

Detector	Reporter	Color	Master Mix	Primer	Probe	Supplies	Batch Number
F1	FAM						
F2	HEX						

Reference Dye

10. Crosstalk Gain Parameter Measurement

The **Crosstalk Gain Parameter Measurement** tool is used to measure crosstalk gain parameters.

Click **Tools ► Measure Crosstalk Gain Parameters ►** the following window will pop up.

The user can add and modify the channels to be tested and dyes according to his needs; upload corresponding reaction plates and operate the experiment. When the experiment is over, the system will automatically save the crosstalk gain parameters.

Experiment Type: **Absolute** Experiment Name: **Crosstalk Gain Measurement**

Experiment Properties

Experiment Name: Comment:

User Name:

Detectors

Detector	Reporter	Color	Master Mix	Primer	Probe	Supplies	Batch Number
F1	FAM						
F2	HEX						

Reference Dye

11. System Maintenance

The **System Maintenance tools** are used for system maintenance.

Click **Tools ► System Maintenance ►** the Password Input box will pop up ► input the correct Password ► conduct the following settings:

- a. Y-axis commissioning
- b. X-axis origin calibration
- c. Machine serial number setting
- d. Photomultiplier setting
- e. Runtime zero clearing
- f. Background measurement
- g. Reference gain measurement
- h. Fluorescence incremental calibration
- i. Firmware Upgrades




Firmware Upgrade tools are used to upgrade the firmware.

Software updates are achieved by connecting to the computer with the RS232 interface supplied with the instrument.

- Set the MODE update switch of the communication box on the back of the instrument to the right hand side ► Update.
- Switch the power on and connect the serial port line.
- The instrument is in update status.
- In the panel, indicator light flashes green and red at the same time, which is normal.
- Click **Tools ► System Maintenance ► Firmware Upgrade ►** the following window will pop up.

The user can:

- a. Select serial ports
- b. Select the BIN file to be upgraded
- c. Upgrade




12. Upgrade Experiment File Format

The **Upgrade Experiment File Format tools** are used to convert old files with the suffix of .fqj or .fqz into new files with the suffix of .fqd.

Click **Tools ► Upgrade Experiment File Format ►** the following window will pop up.


The user can:

- a. Add files to be upgraded
- b. Remove selected files
- c. Select the output directory of new files
- d. Upgrade



13. Ta Calculator

Click **Tools** ► **Ta Calculator** ► the following window will pop up.
Input Forward Primer and Reverse Primer, click Calculate to gain Forward Temperature, Reverse Temperature, Average Temperature and Anneling Temperature.




Chapter 10 Other Functions



1. Instrument Operation

The Instruments operations include **Connect** instrument, **Disconnect** instrument and **Instrument Information**.

1.1 Connect

Click **Instrument** ► **Connect** ► select port number or select automatic port matching.



When the instrument is connected, the icon on the status bar will be ; if the instrument is disconnected, the icon on the status bar will be .

1.2 Disconnect


Click **Instrument** ► **Disconnect** ► disconnect currently connected instrument

1.3 Instrument Information

When the instrument is connected, the user can check the instrument information.

Click **Instrument** ► **Instrument Information** ► the following dialog box will pop up

Instrument information includes instrument serial number, runtime, currently connected ports, and whether an experiment is in operation.



2. Data Query

Data Query is used to query the data already exported to the database.

Click **Data Summary** ► **Data Query** ► the following window will pop up

The user can:

- a. Select database files

- b. Set up query condition
- c. Query
- d. Clear all query conditions

The screenshot shows a software window titled "Data Query". At the top, there is a "Path:" label followed by a text input field and a "Browser..." button. Below this is a "Query Condition" section containing two columns of filters. Each filter consists of a text input field and a dropdown arrow. To the right of these filters are two buttons: "Clear Condition" and "Query".

Below the query condition section is a "Query Result" section. It features a table with the following columns: #, File Name, Sample Id, Sample Name, Test Item, Name, Sex, Age, Case No., Outpatient No., Bed No., Hospital No., Nationality, Sampling Time, Diagnosis, Doctor, Dept., and Test Result. The table body is currently empty.

3. System Help

Click Help ► Help Topics

Chapter 11 Maintenance

1. Regular cleaning

In order to ensure normal operation, detection and use, the instrument needs to be cleaned regularly.

- To clean the outer surface: Clean only with a soft cloth, and if necessary, the cloth may be soaked with alcohol, distilled water or a mild detergent.
- To clean the module wells: Wells may be cleaned with nail wipes which does not bring dust and if necessary, they may be soaked with 95 percent of absolute ethyl alcohol used in medicine or distilled water.

Warning!

1. Before cleaning this instrument, the power supply must be switched off.
2. When cleaning the conical wells of the module, care must be taken to prevent any cleaning agents from dropping into the wells.
3. The surface of this instrument **MUST NOT** be cleaned with corrosive cleaning agents.
4. In order to avoid scratches or damage to the optics in the wells, **NEVER** use sharp or hard objects to clean the wells.

2. Analysis and Troubleshooting


No.	Problem	Possible Cause	Correction
1	The power switch at the rear of the instrument is set ON, but the instrument makes no response.	The run switch at the front of the instrument is not switched on. The internal switching power will have no voltage output.	For working the instrument, its run switch must be pressed and the green indicator lamp should be lit green to indicate energizing of the control system.
2	The display of system parameters menu requires input of "Password".	The system parameters are for instrument manufacturer's internal calibration and require special accession password.	The function is not required for the end user and for calibration contact the manufacturer or supplier's service personnel.
3	After turning the power switch on, the power lamp fails to light up	The run switch is not pressed	This switch is used for temporary turning on/off of the output power and is equivalent to "standby" mode. For prolonged shutdown, the rear power switch should be switched "Off".
		Power is not switched on.	Check and switch on power
		Burnt fuse	Replace fuse (250V 8A Φ 5x20)
		Damaged switch	Replace switch
		Other	Contact the supplier or manufacturer


4	At detecting sample position, the step motor fails to work and the communication fails.	Poor contact or damage of the interface wire	Check, connect or replace interface wire
		The power switch is not turned on or is turned on only after the programme starts running	Turn on the power switch and restart the programme
		The step motor or the drive is damaged	Contact the supplier or manufacturer
		The fixing clamp is not fully inserted	Tighten it into the UNLOCKER port, and switch on power again after shutdown
5	After detecting sample position, the actual temperature displays 0C or 100C	The module temperature sensor is damaged. It accompanies panel red lamp alarm and a software prompt, and the instrument automatically stops running.	Contact the supplier or manufacturer
		The power switch is turned on only after the programme starts running	Switch on the power and restart the programme
		The programme is searching communication port and during this period data would be not sent.	If the trouble still exists after researching, contact the supplier or manufacturer
6	Module temperature heating or cooling rate obviously decreases or temperature control is incorrect.	The ventilation opening is blocked.	Clear the ventilation opening
		Loose connection wire	Contact the supplier or manufacturer
		The refrigerating sheet is damaged	Contact the supplier or manufacturer
		Fan is damaged or fails to run	Contact the supplier or manufacturer
		The temperature sensor is damaged	Contact the supplier or manufacturer
7	The module fails to heat and refrigerate.	The inside of the instrument is damaged	Contact the supplier or manufacturer
		The refrigerating sheet is damaged	Contact the supplier or manufacturer
		During hot-lid heating-up	Waiting until the hot-lid temperature comes to the target value. When stopping running, module temperature holds down 30C automatically.
8	Abnormal temperature or fluorescence curve: straight line or loss of partial data	The running programme is infected by a virus	After removing the virus, re-install the application software
		Computer configuration fails to meet requirements	Configure as per requirements

		or the setup of communication port is not appropriate.	
9	Yellow lamp on panel lights on	The module is not fully pushed in and the optic coupler fails to detect the module	Push in again, if the light is still on, contact the supplier or manufacturer
10	The hot-lid is will not heat	Thermal-sensitive fuse is damaged	Contact the supplier or manufacturer
		Loose plug-pieces	Contact the supplier or manufacturer
		Heating elements of hot-lid is damaged	Contact the supplier or manufacturer
		Temperature sensor of hot-lid is damaged	Contact the supplier or manufacturer
11	Under no test tube state, the fluorescence value difference between wells increases or the background value is very high.	The test tube well or hot-lid is contaminated, or baseline*****.b48 background parameters are set incorrectly.	Eliminate contamination. Each instrument shall correspond to baseline48 document. After perennial use, offset would occur in the optical elements. In this case, contact the manufacturer to re-calibrate the background value.
12	Reagent evaporation	The PCR tube cap does not sealing tightly enough.	Change consumable to one with a tighter fitting cap.
13	Signal crosstalk among channels	Dye signal crosstalk among channels can happen.	You can measure by using “Crosstalk Measurement”, and save parameters to modify.
14	Fluorescence detection value -abnormal	Irradiation by external strong light	Switch off external light source, or remove instrument from external light source
		During a programme run, the hot-lid is opened	Close the hot-lid (detection result unreliable)
		The photo-electric system is damaged	Contact the supplier or manufacturer

Caution: During the warranty period, opening the instrument casing to inspect the internal workings will invalidate the warranty. If any problems should arise please contact the supplier or manufacturer in the first instance.

Appendix: LineGene K series wiring





Note:

A large, empty rectangular box with a thin black border, intended for writing notes. It occupies the central portion of the page.

